



Catalog# ULT286-4A
Ultracold Chest Freezer
4Cu.Ft. (115V)



Catalog# ULT286-7A
Ultracold Chest Freezer
7Cu.Ft. (115V)



Catalog# ULT286-11A
Ultracold Chest Freezer
11Cu.Ft. (115V)



Catalog# ULT286-17A
Ultracold Chest Freezer
17 CuFt. (115V)

Description:

Chest Style Ultra-Low Temperature -86 ° C freezers offer extraordinary value and quality for a variety of applications, these models feature stainless steel interior, EMB condenser and double-seal gasket.

- Designed for optimal efficiency and Temperature stability.
- Dixell Microprocessor controller is simple to use with push-button control and has an easy-to-read operational status display. Large LED show internal temperature, set temperature, ambient temperature with High/Low Temperature alarms
- Patented technology of Autocascade (single compressor) ULT refrigeration, strong cooling, fast freezing and low power consumption!
- Embraco compressor, EBM fan motor from Germany, copper tube evaporator & condenser, stainless steel interior, ensure reliable operation and safety of stored products
- NON-CFC free mixed refrigerants & foaming agent, green and environmentally friendly
- Mounted on Heavy-duty casters for trouble-free/easy movement.

WP SKU	Type	Capacity*	Description	Electrical	Weight	Temp.	Chamber dims (WxDxH)	Exterior dims (WxDxH)
ULT286-4A	Chest	4 cuf (105L)	Ultra-Low Temp. Chest Freezer-86 °C, 115V/60Hz	400W	220 lbs (100 Kg)	-40 to-86 °C	16 x 20 x 22 in. (40.6 x 51 x 56 cm.)	45 x 30 x 37 in. (114.3 x 76.2 x 94 cm.)
ULT286-7A	Chest	7 cuf (280L)	Ultra-Low Temp. Chest Freezer-86 °C, 115V/60Hz	500W	353 lbs (161 Kg)	-40 to-86 °C	28 x 20 x 23 in. (71 x 51 x 58.4 cm.)	55 x 30 x 37 in. (140 x 76.2 x 94 cm.)
ULT286-11A	Chest	11 cuf (308L)	Ultra-Low Temp. Chest Freezer-86 °C, 115V/60Hz	600W	397 lbs (181 Kg)	-40 to-86 °C	40 x 20 x 23 in. (101.6 x 51 x 58.4 cm.)	67 x 30 x 37 in. (170.2 x 76.2 x 94 cm.)
ULT286-17A	Chest	17 cuf (481L)	Ultra-Low Temp. Chest Freezer-86 °C, 115V/60Hz	700W	529 lbs (240.5 Kg)	-40 to-86 °C	51 x 24 x 24 in. (129.5 x 61 x 61 cm.)	78 x 34 x 38 in. (198.1 x 86.4 x 96.5 cm.)

Small Capacity - Chest



Catalog# ULT186-1A
Ultracold Chest Freezer
1Cu.Ft. (115V)



Catalog# ULT186-2A
Ultracold Chest Freezer
2Cu.Ft. (115V)

Description:

Compact Small Volume Chest Style Ultra-Low Temperature -86 ° C freezers offer extraordinary value and flexibility for a variety of applications.

- Designed for optimal efficiency and Temperature stability.
- Dixell Microprocessor controller is simple to use with push-button control and has an easy-to-read operational status display. Large LED show internal temperature, set temperature, ambient temperature with High/Low Temperature alarms
- Patented technology of Autocascade (single compressor) ULT refrigeration, strong cooling, fast freezing and low power consumption!
- Embraco compressor, EBM fan motor from Germany, copper tube evaporator & condenser, stainless steel interior, ensure reliable operation and safety of stored products
- NON-CFC free mixed refrigerants & foaming agent, green and environmentally friendly
- Mounted on Heavy-duty casters for trouble-free/easy movement.

WP SKU	Type	Capacity*	Description	Electrical	Weight	Temp.	Chamber dims (WxDxH)	Exterior dims (WxDxH)
ULT186-1A	Small capacity-Chest	1 cuf (28L)	Ultra-Low Temp. Chest Freezer-86 °C, 115V/60Hz	400W	110 lbs (50 Kg)	-50 to-86 °C	13.8 x 10.6 x 11.2 in. (35 x 27 x 28.5 cm.)	22 x 18.9 x 30.7 in. (56 x 48 x 78 cm.)
ULT186-2A	Small capacity-Chest	2 cuf (50L)	Ultra-Low Temp. Chest Freezer-86 °C, 115V/60Hz	400W	132 lbs (60 Kg)	-50 to-86 °C	17.3 x 13.8 x 12.4 in (44 x 35 x 31.5 cm.)	25.6 x 22 x 34.6 in (65 x 56 x 88 cm.)

Upright Models:



Description:

Upright Style Ultra-Low Temperature-86 °C freezers offer extraordinary value and quality for a variety of applications, these models feature stainless steel interior, EMB condenser and double-seal gasket.

- Designed for optimal efficiency and Temperature stability.
- Dixell Microprocessor controller is simple to use with push-button control and has an easy-to-read operational status display. Large LED show internal temperature, set temperature, ambient temperature with High/Low Temperature alarms
- Patented technology of Autocascade (single compressor) ULT refrigeration, strong cooling, fast freezing and low power consumption!
- Embraco compressor, EBM fan motor from Germany, copper tube evaporator & condenser, stainless steel interior, ensure reliable operation and safety of stored products
- NON-CFC free mixed refrigerants & foaming agent, green and environmentally friendly
- Mounted on Heavy-duty casters for trouble-free/easy movement.



Catalog# ULT386-3A
Ultracold Upright Freezer 3Cu.
Ft. (115V)



Catalog# ULT386-7A
Ultracold Upright Freezer
7 Cu.Ft. (115 Volt)



Catalog# ULT386-25A
Ultracold Upright Freezer
25 Cu.Ft. TWIN COMPRESSORS



Catalog# ULT386-20A
Ultracold Upright Freezer
20 Cu.Ft. (115V)

WP SKU	Type	Capacity*	Description	Electrical	Weight	Temp.	Chamber dims (WxDxH)	Exterior dims (WxDxH)
ULT386-3A	Upright	3 cuf (80L)	Ultra-Low Temp. Chest Freezer-86 °C, 115V/60Hz	400W	176 lbs (80 Kg)	-40 to-80 °C	15 x 16 x 20 in. (38.1 x 40.6 x 51 cm.)	24 x 26 x 50 in. (61 x 66 x 127 cm.)
ULT386-7A	Upright	7 cuf (280L)	Ultra-Low Temp. Chest Freezer-86 °C, 115V/60Hz	400W	220 lbs (100 Kg)	-40 to-80 °C	15 x 16 x 40 in. (38.1 x 40.6 x 102 cm.)	25 x 26 x 71 in. (63.5 x 66 x 180.3 cm.)
ULT386-20A	Upright**	20 cuf (566L)	Ultra-Low Temp. Chest Freezer-86 °C, 115V/60Hz	800W	573 lbs (260.6 Kg)	-40 to-80 °C	23.6 x 27.5 x 51.2 in. (60 x 70 x 130 cm.)	34 x 38 x 78 in. (86.4 x 96.5 x 198.1 cm.)
ULT386-25A	Upright***	25 cuf (708L)	Ultra-Low Temp. Chest Freezer-86 °C, 208/230V/60Hz	Dual Compressor 6.2 Amp	661 lbs (300 Kg)	-50 to-86 °C	35 x 24.4 x 51.2 in. (89 x 62 x 130 cm.)	45.3 x 34.8 x 77.8 in. (115 x 88.4 x 197.6 cm.)

***20 cuf model comes with 4 stainless steel inner doors with 3-shelf configuration

***25 cuf model requires 15.0 A (min.) dedicated grounded circuit. Comes complete with NEMA 6-15P / IEC 10ft. Power cord. Provided with 3 shelves adjustable in 1" increments.